

# OptiClean 900/1300

Vertical Laminar Flow Clean Benches



Heal Force Leads You to Healthier Life



# Enhanced comfort and convenience



5 ° backward-slanted safety grade glass provides more comfortable viewing

Long life HEPA filter offers high separability for particles with high uniformity of air flow

Counterbalanced sliding sash is designed for easy operation with single hand.

Multiple levels of air flow speed are standard configured for wide range of applications

Universal wheels are designed for easy transportation

Advanced door lock guarantees your sample security to prevent unauthorized access

	<b>OptiClean900</b>	<b>Opt</b>
Cleanliness in		ISO Class5, Class 100
Prefilter		High efficiency
Main Supply Filter		HEPA with 99.99% Separability at
Vertical Airflow		0.30 ~ 0.60m/s(Adjustable)
Noise Emission	<65 dB(A)	<65
Fluorescent Light	>300 Lux	>30
Nominal Power	300 VA	550
RMS	<5	<5
External	1050x7	1500
Internal	950 x	1400
Operator No.	One side for one	One
Net Weight	110 Kg	140



General Specifications, OptiClean Vertical Laminar Flow Clean Bench

[www.healforce.com](http://www.healforce.com)

Heal Force Bio-Meditech Holdings Group  
 Heal Force International Trading(Shanghai) Co., Ltd  
 6788 Songze Road, Qingpu District, Shanghai 201706, P.R. China  
 Tel +86 21 62728646 Email [export@healforce.com](mailto:export@healforce.com), Fax +86 21 62728646

